

www.tissyou.com info@tissyou.com

ORTHOPAEDICS

AESTHETIC MEDICINE 🔵

#OP3AB

PMCF Report, 26th February 2024





Andrea Cattalani, Andrea Barbanera

Alessandria, Italy



Hydrolyzed Collagen 5mg/ml - 2ml Medical device based on low molecular weight collagen peptides (LWPs). Ready-to-use solution for the structural strengthening of connective tissues. Produced by Tiss'You.



34 patients 28-69 years



Low Back Pain and facet syndrome



Intra- and peri-articular injections Outpatient treatment with fluoroscopic guidance

0-	-	
	٦J	

Up to 3 months of follow-up

Background

Low back pain (LBP) is a prevalent global symptom, defined as pain in the lower back, extending from the last vertebrae to the buttocks. It is a leading cause of disability, often with nonspecific origins (meaning no identifiable underlying conditions like infection, tumor, or inflammation). Various structures, including intervertebral discs, facet joints, and surrounding tissues, can contribute to LBP. Facet syndrome, characterized by pain originating from the facet joints, is a common source of LBP - particularly in cases of chronic pain. Low molecular weight collagen peptides (LWPs) present a novel bio-

Methods

A cohort of 34 patients, aged between 28 and 69, diagnosed with facet joint syndrome and experiencing lower back pain, participated in this study. The patients were evenly divided into two groups, each undergoing specific treatment protocols. The first group received peri and intra-articular fluoroscopy-guided infiltrations, with 2ml of low molecular weight collagen peptides (LWPs) administered at the initial session (time point 0) and a subselogical approach for LBP treatment. Derived from hydrolyzed bovine collagen, these peptides readily diffuse into surrounding tissues. LWPs can reinforce the extracellular matrix, as they directly support and strengthen weakened connective tissues damaged by degeneration, inflammation, or trauma.

This study aims to evaluate the safety and clinical efficacy of LBP treatment with LWPs compared to cortisone injections, specifically focusing on patients with diagnosed facet arthritis.

quent treatment at 15/20 days. The second group also underwent fluoroscopy-guided procedures, with cortisone injections at time point 0 and a second treatment at 15/20 days. The infiltrations targeted levels such as L4-L5, L5-S1, and the sacroiliac joint. Post-treatment assessments were conducted at 1 and 3 months, utilizing the Visual Analog Scale (VAS) for pain evaluation and the Oswestry Disability Index for functional analysis.

Results

In both study groups, a notable improvement of approximately 1 point on the Visual Analog Scale (VAS) for pain was observed following the initial treatment, administered at time point 0. However, there was no discernible functional enhancement, assessed through the Oswestry Disability Index (ODI), during the second treatment phase at 15/20 days, with no significant differences noted between cortisone and LWPs interventions.

At the 1-month and 3-month follow-ups, a substantial amelioration in both pain levels and functional outcomes was evident compared to pre-treatment conditions in both groups. The improvements were statistically significant when compared to the pre-treatment timeframe, showcasing robust enhancements. Notably, there were no discernible differences between the cortisone and LWPs groups, indicating comparable efficacy in pain relief and functional restoration. The most noteworthy outcomes were observed at the 3-month assessment point, underlining the sustained and progressive nature of the observed improvements.



Figure 1. Mean VAS pain score (**a**) and ODI functional score (**b**) (n=17) before Cortisone and after LWPs treatment. T0 is baseline, TS is secondary treatment at 15/days, T1 is 1-month, T3 is 3-months follow-up.

Errors bars show SEM; #, † p vs. T0 < 0.01

Discussion

The observed comparable efficacy of collagen peptides interventions to cortisone in alleviating pain and enhancing functionality is a promising finding for several reasons. Firstly, LWPs represent a device rather than a pharmaceutical agent, offering a potentially safer alternative for long-term management. Given the well-documented side effects associated with cortisone, the potential of LWPs to achieve similar benefits without these drawbacks is of significant clinical importance.

Moreover, the favorable safety profile of intra and peri-articular injections, including those of LWPs, underscores their tolerability among patients. This suggests a potential advantage in terms of patient acceptance and adherence to treatment regimens.

While our results demonstrate notable improvements within a short timeframe, future investigations with extended follow-up periods would provide valuable insights into the durability and long-term sustainability of the observed effects. The consideration of longer-term outcomes is essential for a comprehensive understanding of the clinical implications of collagen peptide interventions, guiding their potential integration into broader therapeutic strategies for facet joint syndrome and lower back pain.

This report is intended exclusively for medical professionals. The information is provided for reference purposes and should not be considered medical advice. We assume no liability for any consequences arising from its use. The contents are copyrighted by Tiss'You Srl. Any reproduction requires written permission.